

## Procedure for Saliva Analysis

**1.0 Purpose** - This procedure specifies the method of performing saliva analysis in forensic casework.

**2.0 Scope** - This procedure applies to those Forensic Scientists who have been released to do saliva analysis in forensic casework.

### 3.0 Definitions

- **$\alpha$ -Amylase**- an enzyme found in high concentrations in saliva which is responsible for the breakdown of starch into simple sugars.

### 4.0 Equipment, Materials and Reagents

- Sterile disposable scissors or sterile scalpel blade
- Known saliva sample
- RSID kits which contain the test cards and universal buffer
- 1.5 mL centrifuge tube

### 5.0 Procedure

**5.1 Requirements for testing:** Saliva analysis shall be performed on vaginal swabs (including, but not limited to, external vaginal swabs and others similarly identified ), penile swabs, and/or underwear in all cases under the following circumstances:

**5.1.1** When the victim has stated, or is unsure, that cunnilingus or fellatio has occurred.

**5.1.2** Where the alleged suspect has licked, or possibly licked, his/her fingers before penetrating the victim.

**5.1.3** To validate or refute a statement made by a victim, alleged suspect, witness, or other individual involved in the case.

**5.1.4** Saliva analysis will not be performed on rectal swabs or other areas that could be contaminated with feces.

### 5.2 RSID-Saliva

#### 5.2.1 Procedure

**5.2.1.1** Cut a small sample, approximately 0.5 cm<sup>2</sup> (depending on the concentration of the stain), from the evidence sample using sterile disposable scissors or a sterile scalpel blade and place the cutting into a 1.5 mL centrifuge tube.

**5.2.1.2** Add a minimum of 150  $\mu$ L, up to 1 mL, of RSID universal buffer to each sample and mix well. (The amount of buffer added will depend on the sample size; buffer should cover the sample completely.)

**5.2.1.3** Allow the sample to extract for a minimum of 15 min. For weaker or older samples, Forensic Scientists should use a larger quantity of material and/or an extended extraction time to include overnight.

**5.2.1.4** After completing the extraction process, pipette 100 µL of the extracted sample into the sample well on the RSID card.

## **5.2.2 Results**

**5.2.2.1** A positive reaction will have two lines appear in the test window. One line will appear in the area marked "C" for control and one line will appear in the area marked "T" for test. A positive result may be recorded as soon as both of these lines appear, but after no longer than 10 minutes. The lines must be reddish in color.

**5.2.2.2** If a line does not appear in the "T" area within ten minutes, the test is considered negative. A line must appear at the area marked "C" to ensure that the test is working properly.

**5.2.2.3** If no line appears at the area marked "C," the test shall be repeated. If no line is seen in the "C" window in the repeated test, the test is considered inconclusive. If this occurs, the Body Fluid Technical Leader shall be notified as soon as possible. Refer to Forensic Biology Section Administrative Policy and Procedure.

**5.3 Reporting guidelines** – The results statements shall reflect only the work that is performed. Portions of the statements may be omitted to address testing actually performed. This interpretation may include or build upon one (1) or more of the following responses depending on the circumstances of the case and the nature of the examination.

**5.3.1** This phrase shall be used when the RSID Saliva test is negative:

Examination of sample(s) taken from \_\_\_\_ (Item(s) \_\_\_\_), using the RSID Saliva Test, failed to give reactions consistent with the presence of human saliva. No confirmatory saliva testing was performed.

**5.3.2** This phrase shall be used when the RSID Saliva test is positive.

Examination of sample(s) taken from \_\_\_\_ (Item(s) \_\_\_\_), using the RSID Saliva Test, gave reactions consistent with the presence of human saliva. No confirmatory saliva testing was performed.

**5.4 Controls** (to be set up on every case or every batch of cases):

RSID-Saliva: A positive control (applicable body fluid standard), and a negative control (100 µL of universal buffer) shall be run with every case or every batch of cases and the results will be recorded in the case notes as a positive or negative. If a reddish line is seen in the negative control "T" area, the test shall be rerun. If a reddish line appears again in the negative control "T" area, the test shall be considered inconclusive. If this occurs, the Technical Leader shall be notified immediately.

**6.0 Limitations** - Limitations include, but are not limited to, the following: Because the RSID saliva test can cross react with breast milk, the statement used in the report is "gave reactions consistent with the presence of human saliva" instead of "revealed the presence of human saliva."

## 7.0 Safety - N/A

## 8.0 References

Forensic Biology Section Body Fluid training documents

Forensic Biology Section Procedure for Calibration and Maintenance

Forensic Biology Section Procedure for Aseptic Technique and Contamination Control

## 9.0 Records - N/A

## 10.0 Attachments - N/A

| Revision History |                |                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
|------------------|----------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Effective Date   | Version Number | Reason                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| 10/26/2012       | 1              | Original Document - Combined Procedure for Phadebas test and RSID-saliva portion of the RSID Procedure. Added reporting guidelines and allowed for changes to be made by the Forensic Scientist to address the testing actually performed. Removed +/- as a result choice for phadebas test                                                                                                                                                        |
| 12/07/2012       | 2              | 5.2.1 - Added wording to clarify required analysis on an item that has multiple areas being tested.                                                                                                                                                                                                                                                                                                                                                |
| 02/15/2013       | 3              | 5.3.5 – changed “no further confirmatory” to “no confirmatory”<br>6.0 – clarified wording                                                                                                                                                                                                                                                                                                                                                          |
| 09/13/2013       | 4              | 4.0 – removed phadebas equipment and reagents; 5.1 – removed phadebas requirements and reworded what sample types would be tested for saliva; 5.1.1.5 – removed limited sample statement; 5.1.2 – removed phadebas procedure; 5.2.1 – removed RSID-saliva requirements (now in 5.1); 5.3.1, 5.3.2, 5.3.5 – removed reporting guidelines for phadebas testing results; 5.4.1 – removed limitations for phadebas; 6.0- removed controls for phadebas |
|                  |                |                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
|                  |                |                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
|                  |                |                                                                                                                                                                                                                                                                                                                                                                                                                                                    |