



Body Fluid Identification Procedures Index

(Listed by Section)

- 1. Phenolphthalein Procedures**
- 2. Anti-Human Hemoglobin Testing**
- 3. Takayama Procedure**
- 4. Species Origin Determination**
- 5. Ouchterlony Procedure**
- 6. Acid Phosphatase Test**
- 7. Sperm Identification: Christmas Tree Stain**
- 8. ABA Card p30 Test**
- 9. Amylase Test**
- 10. Luminol Test**
- 11. CrimeScope**



1. PHENOLPHTHALIN PROCEDURES

1. Reagents:

Preparation of Stock Solution

Phenolphthalein	4 grams
Sodium Hydroxide pellets	40 grams
Zinc dust	20 grams
Distilled water	1000 ml
Absolute Ethyl Alcohol	Bring up to 1200 ml

Add each reagent of the stock solution to a 5000 ml round bottom refluxing flask. Attach the condensing column to the flask and turn on cold water to column. Heat the flask with an electric heating mantle. Reflux the solution for approximately three hours, until the solution is colorless. After allowing the solution to cool down; decant the liquid into a measured container and use absolute ethyl alcohol to bring the total volume to 1200 ml. Add enough zinc dust to cover the bottom of a dark bottle and pour the phenolphthalein solution into the bottle. Label, date the bottle and store it in the refrigerator at 4°C. Phenolphthalin solution shelf life is 6 months.

An aliquot of phenolphthalin solution is kept at each analyst's bench. A fresh aliquot is prepared the first working day of each month.

Additional reagents needed for the test include:

Absolute Ethanol and 3% Hydrogen Peroxide {prepared from stock 30% H₂O₂}.

2. Sample description:

The blood stain may or may not be visible to the naked eye.

3. Standards and controls:

Standards should include a known blood stain (positive control) and a known blood-free sample (negative control). These controls will be run prior to analysis



and recorded in the laboratory notes.

4. Procedure:

To conduct this test, either rub the suspected stain with a folded piece of filter paper or a clean cotton swab. Add the following reagents in order; one drop of ethanol, one drop of phenolphthalein, and one drop of 3% H₂O₂ onto the sample rubbing. A positive reaction is indicated by the development of a pink color within 5 seconds. Reactions occurring after 5 seconds, or before the addition of the hydrogen peroxide are inconclusive.

NOTE - Phenolphthalein is only a presumptive test for blood and can give reactions for substances other than blood.

5. References:

Blake ET, Dillon DJ (1973) Microorganisms and the presumptive tests for blood. J Police Science and Administration 1: 395-400

Higaki RS, Philp WMS (1976) A study of the sensitivity, stability and specificity of phenolphthalin as an indicator test for blood. Can Soc Forensic Sciences 9: 97-102

Gaensslen RE (1983) Sourcebook in forensic serology, immunology and biochemistry. The National Institute of Justice, Washington, DC p 103-105

6. Safety Precautions

Use care during refluxing of phenolphthalein solution.



2. Anti-Human Hemoglobin Testing

1. Pipette 3 drops (150 μ l) of Hema Trace Extraction Buffer into 1.5ml centrifuge tube.
2. Using sterile scissors, cut at least a 1/4in. fiber from your sample and place into the centrifuge tube using sterile forceps.
3. Then allow the sample to extract for a minimum of 5 min.
4. For weak or older samples, analysts may desire to use a larger quantity of material and a longer extraction time.
5. After completing the extraction process, pipette 3 drops (150 μ l) of liquid buffer into the well marked S on the OneStep ABACard.
6. When there is a positive reaction, two lines will appear, one line in the area marked C for control and one in the area marked T for test. If the reaction is negative then only one line appears in the area marked C. In order to determine that the test is negative a full ten minutes must pass after the liquid is added to the card. If no lines appear the test must be repeated.
7. Beware of a high dose hook effect. If a sample gives a negative or extremely weak positive reaction, and your preliminary (phenolphthalein) testing indicates you would expect a strongly positive result, repeat the testing in the following manner. Add 3 drops of Hema Trace Extraction Buffer to the prior extract tube, pipette the fluid up and down a couple of times to mix the fluid thoroughly, and re-run the test. If necessary, one could repeat the dilutions again in this manner, or do a serial dilution of a new extract.
8. Since validation studies have shown that anti-human hemoglobin reactions were obtained from body fluids other than blood (e.g. urine), then in order for one to identify human blood, both the phenolphthalin and anti-human hemoglobin test must be positive.



3. TAKAYAMA TEST

1. Reagents:

Distilled water	
Sodium Hydroxide pellets	(Sigma S-5881)
Glucose	(Sigma S-5000)
Pyridine	(Fisher P-368)

Preparation of stock solutions:

Dissolve 10 grams of the sodium hydroxide pellets 100 ml of distilled water to prepare a 10% sodium hydroxide solution.

To make a saturated glucose solution add glucose to distilled water in a flask with a magnetic stirrer (on a stir plate) until the glucose precipitates out as a solid.

The 10% sodium hydroxide and saturated glucose stock solutions will be made fresh before use.

Preparation of Working solution:

In a clean test tube add:

Saturated glucose solution	5 drops
10% NaOH solution	5 drops
Distilled water	10 drops
Pyridine	6 drops

Note: The amount of the above solution can be increased by doubling the amount of each chemical added.

2. Sample Description:

A visible dried blood stain.

3. Standards and Controls:



A known dried blood stain should be used as the positive control. A piece of clean cotton cloth should be used as a negative control. These controls will be run with each test and recorded in the laboratory notes.

4. Procedure:

1. Remove a small sample of the stain and place it on a clean microscope slide.
2. Place a coverslip over the sample.
3. The working solution of the Takayama reagent is added to the slide at one edge of the coverslip. The reagent will be drawn under the coverslip by capillary action. Add the working solution until all the area under the coverslip is covered.
4. Let the slide sit at room temperature for a short period of time and then read the slide under the microscope at 40x to 125x magnification via a compound light microscope.
5. If no crystal formation has occurred, heat the slide on a hot plate on a low temperature setting. Be careful not to overheat the slide and dry out the working solution under the coverslip.
6. Once again, read the slide microscopically to see if crystals have formed.
7. A positive reaction will be indicated by pinkish-red rhomboid-shaped crystals.

5. References:

Blake ET, Dillon DJ (1973) Microorganisms and the presumptive tests for blood. J of Police Science 1:397

Gaensslen RE (1983) Sourcebook in forensic serology, immunology and biochemistry. The National Institute of Justice, Washington, DC p85-87

Spaulding RP, Cronin WF Technical and legal aspects of forensic serology: a laboratory manual. p 14-15, 22-24



6. Safety Precautions

Sodium hydroxide pellets are caustic and corrosive.



4. SPECIES ORIGIN DETERMINATION-TUBE TEST

1. Reagents:

Anti-human serum from rabbit (Cappel 5101-1382) or respective antiserum

Rabbit whole serum (Cappel 5012-1380) or respective whole serum

Whole human serum diluted 1:1000 (Cappel 5001-1380) or respective species serum

2. Sample Description:

The sample is an eluate of the suspected blood stain. Soak the blood stain in either saline or distilled water. For best results, try and obtain a slight straw color to give a 1:1000 dilution; however, the reaction may still occur even if no color is noted.

3. Standards and Controls:

The unknown sample will be tested against the species normal serum from which the antisera was produced in, to verify that no non-specific precipitation is occurring. When determining human origin for example, the unknown stain should be checked against rabbit whole serum as well as against anti-human serum. This will assure that there are no false reactions between the unknown sample and the rabbit antigens. To check your antisera to see if it's working, test it against normal human serum (or the respective animal) which has been diluted to yield a straw color, and look for the line of precipitation.

A substrate control will also be run against anti-human serum. The substrate control is from an apparently "unstained" area of the substrate and is run to ensure that material which might react with the anti-serum is not present outside the stained area.



4. Procedure:

1. Using a Pasteur pipette, introduce a small portion of antisera into a clean 6 x 50mm test tube.
2. Using a separate clean pipette introduce a small amount of the whole serum (from the animal species that the antisera was made in) into a separate 6 x 50 mm test tube.
3. Using a third Pasteur pipette, add the unknown eluate into both test tubes layering it carefully to avoid mixing. This will form two layers so that the reaction can occur at the interface.
4. Place the tubes in a holder and do not touch them for about 10-20 minutes.
5. Read the test tubes and observe the interface. If a sharp line appears between the unknown and the antisera and no line appears between the whole serum control and the unknown and both the known human and substrate controls are correct, then a positive reaction is recorded. If a reaction is noted in the whole serum control and the unknown tube, then the test must be repeated. If the reactions are still noted in both tubes, an inconclusive result is reported.

5. References:

Spaulding RP, Cronin WF Technical and legal aspects of forensic serology: a laboratory manual. DOJ/FBI

Tamaki Y, Kishida T (1983) A simple method for the preparation of species-specific anti-human IgG serum. Act Crim Japan 49: 9-12



5. SPECIES ORIGIN - OUCHTERLONY (gel diffusion)

1. Reagents:

Normal saline (0.9 g Sodium Chloride + 100 mls of water)
Agarose Type I (Sigma A6013)
Anti-human sera (Cappel 5101-1382) or respective antisera
Normal rabbit serum (Cappel 5012-1380) or respective whole serum
Normal human serum diluted 1:1000 (Cappel 5001-1380) or respective species serum

Preparation of Agarose Plates:

Add 0.5 g of agar to 100 ml of distilled water and bring the solution to a boil 3 times, then cool. Aliquot 10 mls to the 16 x 125 mm test tubes, place parafilm over the test tubes (or test tubes in a ziplock bag) and store at 4°C. When samples need to be tested, melt the agar in the 10 ml tube and then aliquot out 3 ml of agar into 50 x 9 petri dishes or add the 10 ml of agar to a large square petri dish if numerous samples need to be run.

2. Sample Description:

The sample is an eluate of the suspected blood stain. Soak the blood stain in either saline or distilled water. For best results, try and obtain a slight straw color to give a 1:1000 dilution; however, the reaction may still occur even if no color is noted.

3. Standards and Controls:

The unknown sample will be tested against the whole serum of the animal that the antisera was produced in to detect any non-specific reactions that might occur in the "host" serum, which could lead to a false positive reaction as well as against the species anti-serum of interest. A known blood eluate should be used as a positive control.

A blank extract is prepared from an unstained area of the substrate and tested against the species antiserum of interest, to ensure that material which might react with the anti-serum is not present.




4. Procedure:

After the agar has hardened, cut six holes in a circular fashion and one hole in the center with a cork borer and suction out the agar to form a well.

Use clean Pasteur pipettes and aliquot a small amount of unknown sample eluate, proper antisera, and controls to the appropriate wells. Incubate in a moisture box at 37° C for a minimum of 7 hours, or until standards work. Read the results using a direct light source. If a sharp line appears between the unknown and the antisera and no line appears between the whole serum control and the unknown and both the known human and substrate controls are correct, then a positive reaction is recorded.

An example of how this test may be set-up is as follows:



5. References:

Gaensslen RE (1983) Sourcebook in forensic serology, immunology and

NCSBI MOLECULAR GENETICS SECTION
Technical Procedures Manual
Body Fluid Identification - Ouchterlony
Revision 01



biochemistry National Institute of Justice, Washington, DC pp 224-225

Spaulding RP, Cronin WF Technical and legal aspects of forensic serology: a laboratory manual. DOJ/FBI

Tamaki Y, Kishida T (1983) A simple method for the preparation of species-specific anti-human IgG serum. Act Crim Japan 49: 9-12



6. ACID PHOSPHATASE TEST (Walker Test)

1. Reagents:

Alpha-Naphthyl acid phosphate calcium salt (Sigma N-7250)
Fast Blue Salt B (Sigma D-3502)
Sodium acetate (Fisher S-209)
Distilled water

Preparation of stock solution

Sodium Acetate Buffer:

Add 40g of sodium acetate (CH_3COONa) to two liters of distilled water. Stir the mixture with a magnetic stirrer. Use acetic acid to adjust the pH to 5.0 and store the mixture in a two liter container at 4° C. Shelf life is 2 months.

An aliquot of Sodium Acetate Buffer is kept at each analyst's bench. A fresh aliquot is prepared the first working day of each month.

Fast Blue Working Solution:

Place approximately 5 ml of distilled water in a 13x100 mm tube. Place a small amount of the Fast Blue salt in the water until an opaque yellow color develops.

2. Sample Description

The sample is a cutting of the suspected stain.

3. Standards and Controls

A known seminal stain is used as a positive control. A reagent control is set up where the buffer, alpha-naphthyl acid phosphate calcium salt, and Fast Blue Salt B are used to ensure that the reagents are not giving a false positive result. A substrate control is set up using a control cutting of the same material that the suspected stain is cut from to ensure that something in the material is not causing the test to show a false positive result.

NCSBI MOLECULAR GENETICS SECTION
Technical Procedures Manual
Body Fluid Identification - Acid Phosphatase Test
Revision 01



4. Procedure:

Remove a small portion of the suspected stain area. If the item tested is a pair of panties, at least three cuttings from the crotch area are required. Place each sample into a separately labeled well on a spot plate. Place the positive control, reagent control and cloth controls in their respective wells. Fill the wells about half full with the sodium acetate buffer. Add a small portion of alpha naphthyl acid phosphate calcium salt to each well (a few grains will be sufficient). Mix each well with a separate wooden applicator stick and rotate the plate between 5-15 minutes. After the plate has rotated, add a drop of fast blue solution to each well and look for a color change. A positive reaction will usually turn purple but lighter shades have been noted.

Grade the color change as follows:

Dark Purple - 4+
Purple/Pink - 2+
Pink - 1+
Pale Pink/Yellow - negative (-)

Remember that the AP test is only a preliminary test for semen and that a seminal stain can give a negative AP.

5. References:

American Jurisprudence Proof of Facts (1962) Identification of seminal fluids. Jurisprudence Publishers 12:319-348.

Kind SS (1964) Methods in Forensic Science. A. S. Curry (editor) Interscience Publishers New York 3:267-287.

Schiff Af (1978) Reliability of the acid phosphatase test for the identification of seminal fluid. J Forensic Sciences 23:833-843.

Sensabaugh GF (1975) Genetic and non-genetic variation of human acid phosphatases. Isozymes: Molecular Structure Academic Press New York pp. 367.



6. Safety Precautions

Fast Blue B Salt is an irritant.

Alpha-Naphthyl Acid Phosphate Calcium Salt is an irritant.



7. SPERM IDENTIFICATION - CHRISTMAS TREE STAIN

1. Reagents:

Methanol or ethanol

Preparation of stock Kernechtrot Stain:

Aluminum Sulfate $\text{Al}_2(\text{SO}_4)_3 \cdot 18\text{H}_2\text{O}$ 150 g (Fisher A613)
Nuclear Fast Red 3 g (Sigma N 8002)
Distilled water 3000 ml

In a 4 liter flask, add the $\text{Al}_2(\text{SO}_4)_3$ and water. Stir with a magnetic stirrer until dissolved. Add the Nuclear Fast Red and stir until it appears well dissolved. (Slight warming will enhance the dissolving of the reagent). Filter the solution using a Whatman #1 filter paper and store in a 4 liter container. Shelf life is 1 year at room temperature.

Preparation of stock Picroindigocarmine Stain:

Picric Acid 40g (Fisher A-253)
Indigo Carmine 10g (Sigma I-8130)
Distilled water (3000 ml)

Using great caution, slowly add the picric acid to the distilled water in a 4 liter beaker or flask. Apply low heat and use a magnetic stirrer to get the acid into solution. Make sure the picric acid doesn't get too warm as it may explode under heat. Once the acid has dissolved, add the indigo carmine and stir.

Once the solution cools, filter it with Whatman #1 filter paper and store it in a separate 4 liter container. Shelf life is 1 year at room temperature.

2. Sample Description:

The sample is an air dried smear left on a microscope slide which may have originated from two sources.

1. The air dried smear may have been prepared by a nurse or



- doctor.
2. May be an extract from a suspected semen stain that has been dried.

3. Standards and Controls:

A known sperm slide should be made after the stains have been freshly prepared to ensure that both stains are working correctly. Results of this QC check will be documented in the QC Manual. The QC officer will re-test the stain monthly to see that the stains are still working properly and record this QC check in the QC Manual as well.

4. Procedure:

1. Flame the slide (or place on heat block) briefly before staining it to ensure that the sample is fixed to the slide.
2. Place the slides on a rack and apply the Kernechtrot stain to the slides. Leave the stain on between 15 and 20 minutes.
3. Wash the stain off after this time with methanol or ethanol.
4. Apply the Picroindigocarmine stain to each slide. Leave this stain on only a brief time (no more than 15 seconds).
5. Wash off the stain with methanol or ethanol. Let the slides air dry.
6. Once dry, apply a small amount (a couple of drops) of permount onto the slide and add a 22x50 mm cover slip over the slide.
7. Observe the slide under the microscope at 150x to 350x and confirm the microscopic characteristics of the sperm head at 400x.
8. Record the results.
 1. Spermatozoa have a clear acrosomal cap, a red head and a green tail.
 2. Spermatozoa may be identified without the presence of a tail; but the clear acrosomal cap must be present and clearly visible.

5. References:



Gaensslen RE (1983) Ch 10 in Sourcebook in forensic serology, immunology and biochemistry. National Institute of Justice, Washington DC pp 149-182

American Jurisprudence Proof of Facts (1962) Identification of seminal fluids
Jurisprudence Publishers 12: 319-348

Kaye S (1947) Identification of seminal stains J Criminal Law and Criminology
of Northwestern University 38: 79-83

6. Safety Precautions

Picric Acid is an explosive.



8. ABA Card p30 Test

1. Pipette 5 drops (250 μ l) of HEPES buffered saline (HBS) into 1.5 ml centrifuge tube
2. Using sterile scissors, cut at least a 0.25 cm² cutting from your sample and place into the centrifuge tube using sterile forceps.
3. Then allow the sample to extract for a minimum of 2 hours at room temperature or overnight at 4° C.
4. For weak or older samples, analysts may desire to use a larger quantity of material.
5. After completing the extraction process, centrifuge the samples for 3 minutes at full speed.
6. After centrifugation, pipette 4 drops (~200 μ l) of liquid buffer into the well marked "S" on the OneStep ABACard.
7. When there is a positive reaction, two lines will appear, one line in the area marked "C" for control and one in the area marked "T" for test. If the reaction is negative then only one line appears in the area marked "C". In order to determine that the test is negative a full ten minutes must pass after the liquid is added to the card. If no lines appear the test must be repeated.
8. Beware of a high dose hook effect. If a sample gives a negative or extremely weak positive reaction, and your preliminary (acid phosphatase testing) testing indicates you would expect a strongly positive p30 result, repeat the testing in the following manner. Add 4 drops of HBS to the prior extract tube, pipette the fluid up and down a couple of times to mix the fluid thoroughly, and re-run the test. If necessary, one could repeat the dilutions again in this manner, or do a serial dilution of a new extract.



9. IDENTIFICATION OF AMYLASE

1. Reagents:

Phadebas test tablets (Pharmacia Diagnostics 63105)

Distilled water

0.5 N Sodium Hydroxide solution (Sigma S-5880)

Preparation of stock solution

.5 N Sodium Hydroxide Solution (NaOH):

Add 2.0 g of NaOH pellets to 100 ml of distilled water. Stir with a magnetic stirrer until all of the NaOH pellets are in solution. Store at 4° C. Shelf life is 3 months.

2. Sample Description:

The sample is a small cutting from the suspected stain.

3. Standards and Controls:

A known saliva sample and reagent blank should be set up as positive and negative controls respectively. A small cutting from an unstained area should also be run as a substrate control to ensure nothing in the material is causing a false reaction.

4. Procedure:

1. Make a small cutting of the suspected stain along with your controls and place each cutting in a separate 13x100mm labeled test tube.
2. Add one Phadebas tablet to each tube. Then add 4 ml distilled water to each test tube and vortex the mixture for a few seconds or until the materials are well mixed.
3. Incubate at 37°C for thirty minutes.
4. Remove the test tubes and place 1 ml of 0.5N NaOH in each tube.




5. Vortex each tube for approximately five seconds, then centrifuge each tube for five minutes. This stops the reaction.
6. Results are recorded using a grading system of negative thru 4+, with 4+ being the darkest color. Refer to the color chart below. If area is positive and is going to be analyzed for DNA, a sperm search is still required.

Note: The PCR-STR methods used in the DNA Unit are extremely sensitive and detect as little as 500 picograms of DNA. Analysts are reminded that protein markers used to screen for semen (acid phosphatase and P-30) are more easily degraded than sperm cells, can be affected by various disease states, and are extremely water soluble. For this reason, it is not unexpected that occasionally one will find a sample which yields a negative acid phosphatase result, but is positive for sperm cells.

7. If a result is borderline, it should be confirmed using a spectrophotometer (Spectronic 20). The spectrophotometer should be turned on five minutes prior to the readings. After the five minute warm up, set the spectrophotometer to zero using the left knob on the machine. The wave length should be set to 620nm. Use the reagent blank as your solution blank. Pipet the liquid off the top of the tube into a cuvette, being careful not to get any of the blue solid material with it. Set the full scale using the knob on the right. Discard the fluid in the cuvette, rinse the cuvette with distilled water. Then sample the unknown and record the result. Be sure to rinse the cuvette between samplings with distilled water. If the absorbance is 0.1 or higher and blue in color, then a positive reaction is recorded.



PHADEBAS RESULTS CHART



5. References:

Baxter SJ, Rees B (1975) The identification of saliva in stains in forensic casework. Med Sci Law 15: 37-41.



Caska M (1971) A new type of reagent for the detection of molecular varieties of some hydrolytic enzymes: Detection of (alpha) - amylase isoenzymes. Biochem Journal 121: 575-576.

Gaensslen RE (1983) CH 11 in Sourcebook in forensic serology, immunology and biochemistry National Institute of Justice, Washington, DC pp 183-190.

6. Safety:

Sodium hydroxide is caustic.



10. LUMINOL TEST (Albrecht Reaction)

1. REAGENTS:

5 - Amino - 2,3 dihydro - 1,4 phthalazinedione or 3 aminophthalhydrazide
(0.5 g)
Sodium carbonate (Na_2CO_3) 25g
Sodium perborate ($\text{NaBO}_3 \cdot 4\text{H}_2\text{O}$) 3.5g
Water (distilled) 500ml

WORKING REAGENT:

Weigh out each chemical separately. Place the luminol (5-amino-2,3 dihydro - 1,4 phthalazinedione) and sodium carbonate (Na_2CO_3) in a labeled ziplock bag. Place the sodium carbonate perborate (NaBO_3) in a separate labeled ziplock bag.

Just prior to use, add the contents of both bags into the distilled water and shake until all of the reagents have dissolved and are in solution. Transfer the solution into a hand pump spray bottle.

2. SAMPLE DESCRIPTION:

The sample size may vary from trace amounts which cannot be seen all the way up to large quantities of blood which could be present as well as smears, wipes and other patterns which may be left.

3. STANDARDS AND CONTROLS:

A penny or blood stain is sprayed to ensure that the chemicals are working properly. Note: Luminol is only a presumptive test and can give a reaction for things other than blood.

4. PROCEDURE:

1. Before proceeding make sure permission has been given via consent or a search warrant.
2. Check the spray bottle to ensure a fine mist is being expelled.



3. Usually one starts at a place where the assault has occurred.
Note: Avoid walking over an area that has already been sprayed. Taking this precaution will eliminate unnecessary tracking up of the crime scene.
4. Always spray in front of you and walk backwards while spraying, keeping others behind you. Look for areas where a brightly lit reaction occurs for 5 to 20 seconds.
5. Test these areas by taking a filter paper rubbing of the area and doing a phenolphthalein test on the rubbing.
6. Record only the results that give a positive reaction to both the phenolphthalein test and the luminol test.
7. The luminol reactions obtained may be photographed with a fluorescent ruler present for later comparison.

5. REFERENCES:

Zweidinger RA, Lytle LT, Pitt, CG (1973) Photography of bloodstains visualized by luminol. J. Of Forensic Sciences 18: 296-302.

Proescher F, Moody Am (1939) Detection of blood by means of chemiluminescence. The Journal of Laboratory and Clinical Medicine 1183-1189.

Blake ET, Dillon DJ (1973) Microorganisms and the presumptive tests for blood. J Police Science and Administration 1: 395-400

6. SAFETY PRECAUTIONS

Use gloves when handling powder and liquid solutions.



11. CRIMESCOPE

1. Suggested Uses: This procedure is used to examine items of evidence for the presence of biological evidence which may not be visible to the naked eye.
2. Operating Instructions:
 1. Start up procedures:
 1. Set the wheel to “ White light “ as a default when not using the unit.
 2. Turn on the main switch (in the back). Check with your hands that there is air circulation on both exhaust top holes and on the fan located above the light guide connection.
 3. Aim the wheel towards a wall and turn on the lamp switch (on the front of the unit). Within 1 - 2 minutes you should see a bright spot. Do not look at the spot for too long as it is extremely intense. If the spot is weak, verify that the intensity knob on the front of the unit (located below the light guide connection) is open all the way.
 4. Scan your filters for the different colors available. Do not look into the light source.
 5. Start using protective goggles after checking colors.
 6. Turn out lights and get room as dark as possible
 2. Examination procedures
 1. Scan the item of evidence with the light beam.
 2. Collect and/or make notes as appropriate.
 3. Refer to the information below for wave

NCSBI MOLECULAR GENETICS SECTION**Technical Procedures Manual****Body Fluid Identification - Operation of the MCS - 400 Mini - Crimescope****Revision 01**

length and goggle use:

BODY FLUID DETECTION
Semen, Saliva, Urine, Sweat, Vaginal fluids, Feces...
Blood does not fluoresce: it absorbs at 415 nm and reflects at 254 nm.
For dark surfaces and for saliva use "UV" and Clear/Yellow goggles.
On most clothes and rugs use "445/455/CSS/515" and Orange goggles. Use 515 on white clothes .
Refer to operators manual for further applications.

3. Shut down procedures:

1. Turn off lamp (front switch) and let the fans run for 3 - 5 minutes.
2. Turn off main switch (back switch).
3. Do not restart the lamp when the bulb is still hot. An arcing noise may be heard when re-starting too early.
4. If power is lost in the building, turn the lamp off but leave the fan switch on so the lamp can continue to cool when the power is restored.

4. Safety Precautions:

Always wear protective goggles when operating the Mini-Crimescope.

NCSBI MOLECULAR GENETICS SECTION**Technical Procedures Manual****Body Fluid Identification - Operation of the MCS - 400 Mini - Crimescope****Revision 01**

Revision History		
Effective Date	Revision Number	Reason
	00	Original Documents
November 22, 2001	01	1) Procedure Updates 2) Collation of all Body Fluid Procedures into one protocol