

DRUG CHEMISTRY SECTION TECHNICAL PROCEDURE MANUAL		
Procedure C-03	Thin Layer Chromatography 4:1 Developing Solvent	
Effective Date:	November 20, 2006	Page 1 of 3

**Name of Procedure:**

Thin-Layer Chromatography  
4:1 Developing Solvent

**Suggested Uses:**

The 4:1 Thin-Layer Chromatography Developing Solvent is used to separate and identify organic compounds. Suggested uses of this TLC solvent are the cannabinoids found in marijuana.

**Apparatus Needed to Perform Procedure Including Preparation of Reagent:**

Thin-Layer Chromatography Plates  
     Silica gel GF - fluorescent indicator  
     Thickness: 250 microns  
 Developing tank  
 Micro pipets  
 UV light (long and short wave)  
 Fume hood  
 Gloves  
 Eye protection  
 Laboratory coat  
 Graduated cylinders  
 Reagent bottle  
 Funnel  
 Porcelain spot plate  
 Spatula  
 Ethyl ether  
 Hexane

**Formula for Preparing Reagent:**

1. Measure out 40 milliliters of hexane and 10 milliliters of ethyl ether.
2. Pour measured solvents into reagent bottle and mix well.
3. Properly label reagent bottle.

DRUG CHEMISTRY SECTION TECHNICAL PROCEDURE MANUAL		
Procedure C-03	Thin Layer Chromatography 4:1 Developing Solvent	
Effective Date:	November 20, 2006	Page 2 of 3

### **Quality Control Check:**

A quality control check of this reagent will be performed using a known standard of hash and following the application procedure listed below.

### **Expiration Date of Chemical Reagent:**

The solvent can be used until depletion provided it is stored in an airtight reagent bottle in a cool place.

### **Application of Procedure on Evidence:**

1. A sample of the unknown (approximately 1-2 milligrams) is placed in the well of a porcelain spot plate and several drops of a suitable solvent are added to dissolve sample.
2. A known standard (approximately 1-2 milligrams) is also dissolved in the well of a spot plate with several drops of a suitable solvent. (Known standard is selected after visual observation and/or a series of preliminary tests).
3. With a capillary pipet, several microliters of the unknown and standards solutions are placed side by side approximately 5 millimeters from the bottom of the TLC plate.
4. Add 4:1 solvent to the developing tank to a depth of approximately 2 millimeters and allow several minutes for atmosphere to equilibrate.
5. Allow all the spotting solvent to evaporate from the TLC plate.
6. Place the TLC plate in the developing tank and close the lid.
7. Allow the TLC plate to develop to the top of the plate.
8. Remove the TLC plate from the developing tank and allow the solvent to dry.
9. Visualize with UV light or an appropriate visualizing reagent.
10. Record results.

DRUG CHEMISTRY SECTION TECHNICAL PROCEDURE MANUAL		
Procedure C-03	Thin Layer Chromatography 4:1 Developing Solvent	
Effective Date:	November 20, 2006	Page 3 of 3

### **Safety Concerns:**

Always wear eye protection, gloves and a laboratory coat when preparing this TLC solvent.

### **Literature References:**

Randerath, Kurt, **Thin-Layer Chromatography**, New York, Academic Press, 1968

Moffat, A.C., **Clarke's Isolation and Identification of Drugs**, 2nd Ed., The Pharmaceutical Press, 1986, pp. 166-177.

Parker, James M. and Fiske, Harold L., ■Thin Layer Chromatography of Marihuana,■  
**Journal of Forensic Sciences**, Vol. 24, No. 4, 1979, pp. 842-846.