DRUG CHEMISTRY SECTION TECHNICAL PROCEDURE MANUAL				
Procedure D-07	Extraction and Separations			
Separation of Organic Acids and Bases by Solvent Wash				
Effective Date:	November 20, 2006	Page 1 of 2		

Name of Procedure:

Extractions and Separations Separation of Organic Acids and Bases by Solvent Wash

Suggested Uses:

Organic acids and bases are commonly encountered in mixtures with sugars (mannitol, inositol) or other common diluents. The differences in solubility between the organic acids and bases and these diluent materials can be used to separate the organic components for infrared analysis.

Apparatus Needed to Perform Procedure Including Preparation of Reagents:

Fume hood Eye protection Gloves Laboratory coat Small beaker Filter paper Chloroform Ethyl Ether Hexane Methylene Chloride Acetone Spatula, small Heat source

Application of Procedure on Evidence:

- 1. Place 10-30 milligrams of sample in filter paper over small beaker.
- 2. Wash sample with several small portions of suitable solvent.
- 3. Evaporate solvent over heat source in a fume hood to yield compounds.

Safety Concerns:

Ethyl ether, hexane and acetone are extremely flammable. Chloroform and methylene chloride should be used in a well ventilated area or in a fume hood.

This document is not controlled if printed.

DRUG CHEMISTRY SECTION TECHNICAL PROCEDURE MANUAL				
Procedure D-07	Extraction and Separations			
Separation of Organic Acids and Bases by Solvent Wash				
Effective Date:	November 20, 2006	Page 2 of 2		

Other:

This procedure can also be used to remove controlled substances from commercial preparations.

Examples: Diazepam can be removed from tablets with an acetone or ethyl ether wash of the crushed tablet and methylphenidate can be removed using chloroform.

Literature References:

Forensic and Analytical Chemistry of Clandestine Phenethylamines, CND Analytical, Inc., 1994, pp. 6-7.

Moffat, A. C. Ed., <u>Clarke s Isolation and Identification of Drugs</u>, 2nd Ed., The Pharmaceutical Press, London 1986.