	DV OF OTION TEOLINION DECOFERING MAN	11.1.4.1	
DRUG CHEMISTRY SECTION TECHNICAL PROCEDURE MANUAL			
Procedure D-16	Extraction and Separations		
	Extraction of Lysergic Acid Diethylamide (LS	SD)	
Effective Date:	November 20, 2006	Page 1 of 2	

Name of Procedure:

Extractions and Separations
Extraction of Lysergic Acid Diethylamide (LSD)

Suggested Uses:

This procedure is used to extract LSD from capsules, liquid, tablets or blotter paper.

Apparatus Needed to Perform Procedure Including Preparation of Reagent:

Fume hood
Gloves
Eye protection
Laboratory coat
Sodium Bicarbonate
Chloroform
pH Test paper
Sodium Sulfate, anhydrous
Pipets, glass, disposable
Pipet bulb
UV Light source

Application of Procedure on Evidence:

- 1. Dissolve sample in water.
- 2. Add sodium bicarbonate until pH is 8. Check with test paper.
- 3. Extract with chloroform and evaporate over heat source for further analysis.

Safety Concerns:

Chloroform should be used in a well ventilated area or in a fume hood. LSD should not come into contact with skin.

DRUG CHEMISTRY SECTION TECHNICAL PROCEDURE MANUAL		
Procedure D-16	Extraction and Separations	
	Extraction of Lysergic Acid Diethylamide (LS	SD)
Effective Date:	November 20, 2006	Page 2 of 2

Literature References:

Zacharias, D., ■Infrared Identification of LSD in Capsules, Liquid, Sugar Cubes and Tablets■, MICROGRAM, Vol. I, No. 4, Jan. 1968, pp. 23-24.