

Drug Chemistry Section
Drug Chemistry Procedure Manual
Effective Date: March 15, 2003

Modification of J-6
Prepared by R. W. Waggoner, Jr.
Approved by D. J. Koontz
Supercedes: October 1, 2001

Name of Procedure:

Toxicology
THC and THC-COOH Extraction Procedure Using United Chemical Technologies Clean Screen Extraction Columns[®]

Suggested Uses:

This procedure is an extraction of delta-9-tetrahydrocannabinol (THC) and 11-nor-delta-9-tetrahydrocannabinol-9-carboxylic acid (THC-COOH) from blood using United Chemical Technologies Clean Screen Extraction Columns[®]. The procedure is designed to extract THC and THC-COOH for confirmation by mass spectrometry. The procedure targets THC and THC-COOH in a 1.0 mL blood sample. Calibration standards may be utilized for quantitative analysis.

Items Used to Perform Procedure:

Test tubes, 16 x 125, 13 x 100, 12 x 75
Test tube caps or stoppers
Vortexer
Centrifuge
pH paper
pH meter
Eppendorf pipettes
Pipet tips
Sonicater
Volumetric flasks
World Wide Monitoring Clean Screen Extraction Columns[®] CSTHC203
Zymark RapidTrace SPE Workstation
Zymark TurboVap LV
Nitrogen at 30 psi

Reagents Used:

Acquired drug standards
 Delta-9-tetrahydrocannabinol (THC)
 11-nor-delta-9-tetrahydrocannabinol-9-carboxylic acid (THC-COOH)

Drug Chemistry Section
Drug Chemistry Procedure Manual
Effective Date: March 15, 2003

Modification of J-6
Prepared by R. W. Waggoner, Jr.
Approved by D. J. Koontz
Supercedes: October 1, 2001

Acquired deuterated drug standards

Delta-9-tetrahydrocannabinol (THC)-D₃

11-nor-delta-9-tetrahydrocannabinol-9-carboxylic acid-D₃

Deionized water

Methanol

Hexane

Ethyl Acetate

Acetonitrile

BSTFA with 1% TMCS (N,O-bis(trimethylsilyl)trifluoroacetamine with 1% trimethylchlorosilane)

Hydrochloric acid, 0.1M

- Add 4.2 mL concentrated hydrochloric acid to approximately 400 milliliters of deionized water.
- Dilute to 500 mL with deionized water. Mix.
- Label with expiration date. Stability: 6 months at room temp.

Sodium Acetate buffer, 0.1M, pH 4.5

- Dissolve 2.93 g of Sodium acetate trihydrate in 400 mL deionized water.
- Add 1.62 mL of glacial Acetic Acid.
- Dilute to 500 mL with deionized water.
- Adjust pH to 4.5 +/- 0.1 with 100 mM sodium acetate or 100 mM acetic acid.
- Label with expiration date. Stability: 6 months at room temp.

Internal Standard Solution

- Internal standard: Prepare an internal standard solution that contains 1.0 microgram/mL (1,000 ng/mL) of THC-D₃ and THC-COOH-D₃. Label the flask "Cannabinoid Internal Standard". Include on the label initials, date prepared, and the expiration date which is determined from earliest date listed on the acquired drug standards.
Example: From separate drug standards containing 100 micrograms/mL of THC-D₃ and THC-COOH-D₃, transfer 1.0 mL from each standard to the same 100 mL volumetric flask. Dilute the flask to volume with methanol.

Calibration and Verification Solutions

- Cannabinoid Calibrator Solution: from the individual drug standards prepare a solution that contains 1.0 microgram/mL (1,000 ng/mL) of THC and THC-COOH. Label this solution "Cannabinoid Calibrator Solution" and include initials, the

date prepared and the expiration date which is determined from the earliest date listed on the acquired drug standards.

Example: Dilute 1.0 mL of a 1.0 mg/mL THC standard to 10 mL with methanol. Place 1.0 mL of this solution in a 100 mL volumetric flask and add 1.0 mL of a 100 microgram/mL THC-COOH standard. Dilute the flask to volume with methanol.

- b. Cannabinoid Verification Solution must be prepared using the same procedure to prepare the calibrator solution.
- c. 50 ng/mL cannabinoid calibration/verification standard: Add 50 microliters of the cannabinoid calibration/verification solution to 1.0 mL of drug free blood. Add 100 microliters of the internal standard solution. Vortex the mixture.
- d. 100 ng/mL cannabinoid calibration/verification standard: Add 100 microliters of the cannabinoid calibration/verification solution to 1.0 mL of drug free blood. Add 100 microliters of the internal standard solution. Vortex the mixture.
- e. 200 ng/mL cannabinoid calibration/verification standard: Add 200 microliters of the cannabinoid calibration/verification standard to 1.0 mL of drug free blood. Add 100 microliters of the internal standard solution. Vortex the mixture.

Procedure:

1. Blood Sample Preparation:
 - a. To 1 mL of blood add 1 mL of acetonitrile and 100 microliters of the internal standard.
 - b. Mix/Vortex and let stand for 5 minutes.
 - c. Sonicate sample for at least 1 minute. (This step is not required)
 - d. Mix/Vortex samples.
 - e. Centrifuge for 10 minutes at >2000 RPM
 - f. Decant liquid portion of the sample into 5 mL of 0.1M Sodium Acetate buffer pH 4.5. Mix.
2. Extraction Procedure
 - a. Condition column with 3 mL MeOH.
 - b. Condition column with 3 mL water.
 - c. Condition column with 1 mL 0.1M HCl.
 - d. Load sample onto column.
 - e. Rinse column with 2 mL water.
 - f. Rinse column with 2 mL (0.1M HCL : Acetonitrile (70:30))

- g. Dry column for 3 minutes with nitrogen
- h. Rinse column with 0.2 mL Hexane.
- i. Collect cannabinoids with 3 mL of (Hexane : Ethyl acetate (75:25))

Reagent flow must be between 1 to 15 mL per minute, except for elution reagent and samples which cannot be more than 5 mL per minute.

Post Extraction Procedure:

1. Evaporate the solvent from the collection test tube.
2. If desired, prior to derivatization the eluate may be transferred from the collection test tube to a vial using a hexane : ethyl acetate (75:25) mixture.
3. Derivatize by adding 50 microliters BSTFA with 1% TMCS to the vial or collection test tube and capping. Mix and heat the vial or collection test tube at 80° C for 30 minutes. Remove from heat source and allow the vial or collection test tube to cool before analysis.
4. Examine the extract utilizing gas or liquid chromatography/mass spectrometry.

Quality Control:

Quality control is verified for each extraction by utilizing the internal standard. For each set of extractions a blood blank must be extracted as a negative control. At least one verifier must be extracted as a positive control. Reported quantitative values must be based on a calibration curve using at least three calibration standards, and at least one verification standard must be extracted to verify the calibration curve. The verification standard must quantify within +/- 20% of the target for quantitative values to be reported.

Personnel Safety Concerns:

When working with biohazardous samples use protective measures, such as gloves, eye protection, and work with the samples in a biosafety hood.

BSTFA with 1% TMCS should be handled in a fume hood, with gloves, and eye protection.

Maintenance:

Zymark: RapidTrace SPE Workstation

- a. Check reagent levels daily before using.
- b. Clean protein build-up when needed by passing about 7.5 mL 2N NaOH and 2N HNO₃ through the cannula to the aqueous waste.

Comments:

For an explanation of the operation of the RapidTrace refer to the operation manual.

Literature References:

RapidTrace SPE Workstation Installation and Quick Reference Manual, revision 0, Zymark Co., 1995.

Clean Screen® Extraction Column Applications Manual, United Chemical Technologies, Inc., Bristol, PA.