

TABLE OF CONTENTS

Capillary Electrophoresis Information	2
Part 1: Running the CE System	3
1.1 CE System Pre-operation Procedure.....	3
1.2 CE System Parameters and Power Up	3
1.3 PROGRAM 6 Parameters.....	4
1.4 SCHEDULE Parameters	5
1.5 Automatic Injection Using CE System and AI-450 Software.....	6
1.6 Manual Injection Using CE System.....	6
Part 2: Running the AI-450 Chromatography Software	7
2.1 Method Program.....	7
2.2 Schedule Program.....	7
2.3 Run Program	7
Part 3: Preparation of 10 ppm Individual Cation Standards.....	8
3.1 Ammonium Standard from Ammonium Nitrate (NH_4NO_3)	8
3.2 Barium Standard from Barium Chloride (BaCl_2)	8
3.3 Potassium Standard from Potassium Sulfate (K_2SO_4)	8
3.4 Magnesium Standard from Magnesium Nitrate ($\text{Mg}(\text{NO}_3)_2$).....	9
3.5 Sodium Standard from Sodium Chlorate (NaClO_2).....	9
3.6 Lead Standard from Lead Nitrate ($\text{Pb}(\text{NO}_3)_2$).....	9
3.7 Cation Standard Preparation	9
Part 4: Conditioning Capillary	10
4.1 Pressure Injection of 0.1 N NaOH.....	10
4.2 Rinsing Capillary with dH_2O	11
Part 5: Cation Analysis of Unknown Sample	12
5.1 Procedure for Cation Analysis of Unknown Sample	12
Part 6: Storing Capillary	12
6.1 Rinsing Capillary with dH_2O	12
6.2 Drying Capillary with Air.....	12

Capillary Electrophoresis Information

Dionex Capillary Electrophoresis System

Serial No. 941206

System Unit No. 16451

PO No. 0916000724

AI-450 Chromatography Software

Release 3.32

Dionex Advance Computer Interface

Model ACI-1

Serial No. 941408

NEC Technologies Monitor

MultiSync XV₁₅₊

Serial No. 6505994EF

Model No. JC-1571VMA

Dionex IonPhor Cation DDP Electrolyte Buffer

5 mM Dimethyldiphenylphosphonium
Hydroxide (DDP)

2 mM 18-Crown-6

6 mM 2-Hydroxyisobutyric Acid

pH 4.4 - 4.8

Product No. 46071

Dell Optiplex Computer

Model 466/L

Hewlett Packard Desk Jet PLUS Printer

Part 1: Running the CE System

1.1 CE System Pre-operation Procedure

- 1.11 Check buffer bottles to ensure enough buffer is present.
- 1.12 Buffer A is Millipore dH₂O.
- 1.13 Buffer B is Dionex IonPhor Anion PMA Electrolyte Buffer.
- 1.14 Buffer C is Dionex IonPhor Cation DDP Electrolyte Buffer.
- 1.15 Check waste vial in front of Buffer C, and empty if more than half full.
- 1.16 Proceed to 1.2.

1.2 CE System Parameters and Power Up

- 1.21 Turn on instrument in the following order, helium, CE System, ACI, computer, monitor, printer.
- 1.22 Set WAVELENGTH nm to 215.
- 1.23 The top row of parameters are set to LOCAL, ABSORB, UV LAMP to START, and VIS LAMP to OFF.
- 1.24 The OUTPUT RANGE is set to 0.01.
- 1.25 Press ENTER after each selection when using number pad or editing a program/schedule or the option will not be changed.
- 1.26 Proceed to 1.3.

1.3 PROGRAM 6 Parameters

1.31 Select 1 from the MAIN MENU to EDIT (view) PROGRAM 6.

1.32 Select 5 to Edit/Review PROGRAM NUMBER 6.

1.33 Check the following parameters set for cation analysis:

Polarity	(+); detector side (-)		
Air cooling	off		
Rinse Destination	6 seconds	Buffer C	
Rinse Capillary	120 seconds	Destination Buffer	
Rinse fixed Source	5 seconds	Buffer C	
Injection Method	gravity, 100mm for 30 seconds		
Control Mode	voltage		
Current Limit	200 μ A		
Power Limit	2000 mW		
0.0 min	0 V	Relay On: 1	Source Vial Fixed
0.1 min	2000 V	Relay Off: 1	Source Vial Fixed
1.0 min	20000 V		Source Vial Fixed
End of Program	9.0 min		

1.34 Push ESCAPE to exit EDIT mode.

1.35 Select 1 of 3 options, SAVE, SAVE AS NEW NUMBER, or DISCARD.

1.36 Check to make sure voltage wand is in (+) position.

1.36 Proceed to either 1.4 or 1.6.

1.4 SCHEDULE Parameters

- 1.41 Select 3 from the MAIN MENU to Edit/Review a SCHEDULE.
- 1.42 Select 3 to Edit/Review SCHEDULE NUMBER 3.
- 1.43 Enter the number of ITERATIONS (injections) for each sample vial for STEP 1.
- 1.44 Enter the number of the vials to be injected for STEP 1.
- 1.45 Enter the number of the PROGRAM to be used for STEP 1.
- 1.46 Repeat for STEPS 2...n if needed.
- 1.47 Push ESCAPE to exit SCHEDULE mode.
- 1.48 Select 1 of 3 options, SAVE, SAVE AS NEW NUMBER, or DISCARD.
- 1.49 Proceed to 2.1.

1.5 Automatic Injection Using CE System and AI-450 Software

1.51 SCHEDULE must be “loaded” on the software and ready to record data before proceeding.

1.52 Select 6 from the MAIN MENU to run a SCHEDULE in connection with the AI-450 Software.

1.53 Enter SCHEDULE number. Press ENTER.

1.54 Select 1 to START and RUN ENTIRE SCHEDULE NOW.

1.6 Manual Injection Using CE System

1.61 Select 5 from the MAIN MENU to run a PROGRAM.

1.62 Select desired PROGRAM NUMBER.

1.63 Select either RUN/RINSE/INJECT/EVENTS/next vial.

Part 2: Running the AI-450 Chromatography Software

2.1 METHOD Program

- 2.11 Select METHOD from the MAIN MENU.
- 2.12 All methods beginning with CZE215. . . are valid methods for cation analysis.
Example: CZE2159.met --- This method is set for a 9 minute run. The cations validation requires a 9 minute run.
- 2.13 Proceed to 2.2

2.2 SCHEDULE Program

- 2.21 Select SCHEDULE from the MAIN MENU.
- 2.22 Enter the sample name, method, and data file name for each injection.
Note: a method chosen must already be saved or the computer will not comprehend and will not save.
- 2.23 Select Save As. . . under the File menu.
- 2.24 Save the SCHEDULE according to the date plus a letter, i.e. "030399a.sch." Click OK.
- 2.25 Proceed to 2.3.

2.3 RUN Program

- 2.31 Select RUN from the MAIN MENU.
- 2.32 Select SCHEDULE from the Load. . . menu.
- 2.33 Enter the SCHEDULE name or select by scrolling to the schedule and double clicking.
Click OK.
- 2.34 Select "Upon receiving signal at interface". Click OK. Proceed to 1.5.

Part 3: Preparation of 10 ppm Individual Cation Standards

* All cation samples and standards must be prepared using plasticware.

3.1 Preparation of 10 ppm Ammonium Standard from Ammonium Nitrate (NH_4NO_3)

- 3.11 Weigh out 0.02217g NH_4NO_3 and place in 13X100mm disposable plastic test tube.
- 3.12 Add 5ml deionized water (dH_2O).
- 3.13 Cover with Parafilm. Vortex.
- 3.14 Place 100 μl of previous $\text{NH}_4\text{NO}_3(\text{l})$ into a 12x75mm polypropylene tube with cap and add 900 μl dH_2O .
- 3.15 Vortex.
- 3.16 Repeat 3.14 and vortex.

3.2 Preparation of 10 ppm Barium Standard from Barium Chloride (BaCl_2)

- 3.21 Weigh out 0.0314g BaCl_2 and place in 13X100mm disposable plastic test tube.
- 3.22 Add 5ml deionized water (dH_2O).
- 3.23 Cover with Parafilm. Vortex.
- 3.24 Place 100 μl of previous $\text{BaCl}_2(\text{l})$ into a 12x75mm polypropylene tube with cap and add 900 μl dH_2O .
- 3.25 Vortex.
- 3.26 Repeat 3.24 twice, vortexing between each dilution.

3.3 Preparation of 10 ppm Potassium Standard from Potassium Sulfate (K_2SO_4)

- 3.31 Weigh out 0.0111g K_2SO_4 and place in 13X100mm disposable plastic test tube.
- 3.32 Add 5ml deionized water (dH_2O).
- 3.33 Cover with Parafilm. Vortex.
- 3.34 Place 100 μl of previous $\text{K}_2\text{SO}_4(\text{l})$ into a 12x75mm polypropylene tube with cap and add 900 μl dH_2O .
- 3.35 Vortex.
- 3.36 Repeat 3.34 and vortex.

3.4 Preparation of 5 ppm Magnesium Standard from Magnesium Nitrate ($\text{Mg}(\text{NO}_3)_2$)

- 3.41 Weigh out 0.0527g $\text{Mg}(\text{NO}_3)_2$ and place in 13X100mm disposable plastic test tube.
- 3.42 Add 5ml deionized water (dH_2O).
- 3.43 Cover with Parafilm. Vortex.
- 3.44 Place 100 μl of previous $\text{Mg}(\text{NO}_3)_2(\text{l})$ into a 12x75mm polypropylene tube with cap and add 900 μl dH_2O .
- 3.45 Vortex.
- 3.46 Repeat 3.44 and vortex.
- 3.47 Take last dilution and add 1000 μl dH_2O to prepare 2x dilution and vortex.

3.5 Preparation of 10 ppm Sodium Standard from Sodium Chlorate (NaClO_3)

- 3.51 Weigh out 0.02315g NaNO_3 and place in 13X100mm disposable plastic test tube.
- 3.52 Add 5ml deionized water (dH_2O).
- 3.53 Cover with Parafilm. Vortex.
- 3.54 Place 100 μl of previous $\text{NaNO}_3(\text{l})$ into a 12x75mm polypropylene tube with cap and add 900 μl dH_2O .
- 3.55 Vortex.
- 3.56 Repeat 3.54 twice, vortexing between each dilution.

3.6 Preparation of 10 ppm Lead Standard from Lead Nitrate ($\text{Pb}(\text{NO}_3)_2$)

- 3.61 Weigh out 0.0799g $\text{Pb}(\text{NO}_3)_2$ and place in 13X100mm disposable plastic test tube.
- 3.62 Add 5ml deionized water (dH_2O).
- 3.63 Cover with Parafilm. Vortex.
- 3.64 Place 100 μl of previous $\text{Pb}(\text{NO}_3)_2(\text{l})$ into a 12x75mm polypropylene tube with cap and add 900 μl dH_2O .
- 3.65 Vortex.
- 3.66 Repeat 3.64 three times, vortexing between each dilution.

3.7 Preparation of Cation Standard

- 3.71 Place 100 μl of each individual cation standard into a 12x75mm polypropylene tube with cap.
- 3.72 Add an extra 100 μl of Ba^+ and Pb^{2+} , then qs. to 1000 μl .
- 3.73 Vortex.

Part 4: Conditioning Capillary

- * Only 0.1 N NaOH is sufficient to condition capillary for cation method.
Only condition a new capillary prior to initial use.
Capillary does not need to be conditioned after initial installation.

5.1 Pressure Injection of 0.1 N NaOH

5.11 Edit (view) PROGRAM 20 to ensure the following parameters are set:

Polarity	(+); detector side (-)		
Air cooling	off		
Rinse Destination	6 seconds	Buffer C	
Rinse Capillary	120 seconds	Destination Buffer	
Rinse fixed Source	5 seconds	Buffer C	
Injection Method	pressure for 240 seconds		
Control Mode	voltage		
Current Limit	200 μ A		
Power Limit	2000 mW		
0.0 min	0 V	Relay On: 1	Source Vial Fixed
0.1 min	0 V	Relay Off: 1	Source Vial Fixed
1.0 min	0 V		Source Vial Fixed
End of Program	1.0 min		

5.12 Select 5 from the MAIN MENU.

5.13 Select PROGRAM 20.

5.14 Select INJECT from the RUN/RINSE/INJECT/EVENTS/next vial.

5.15 Leave NaOH in capillary for 10 minutes before rinsing.

5.2 Rinsing Capillary with dH₂O

5.21 Edit (view) PROGRAM 3 to ensure the following parameters are set:

Polarity	(+); detector side (-)
Air cooling	off
Rinse Destination	9 seconds Buffer A
Rinse Capillary	120 seconds Destination Buffer
Rinse fixed Source	7 seconds Buffer A
Injection Method	pressure 5 for 5 seconds
Control Mode	voltage
Current Limit	200 μ A
Power Limit	2000 mW
0.0 min	0 V Relay On: 1 Source Vial Fixed
0.1 min	0 V Relay Off: 1 Source Vial Fixed
1.0 min	0 V Source Vial Fixed
End of Program	3.0 min

5.22 Select 5 from the MAIN MENU

5.23 Select PROGRAM 3 to rinse capillary with Buffer A.

5.24 Select RINSE from the RUN/RINSE/INJECT/EVENTS/next vial.

5.25 Repeat rinse a total of 3 times.

Part 5: Cation Analysis of Unknown Sample

- 5.1 Procedure for Cation Analysis of Unknown Sample
 - 5.11 Inject a cation standard according to cation procedure.
 - 5.12 Inject the unknown according to cation procedure.

Part 6: Storing Capillary

- 6.2 Rinsing Capillary with dH₂O
 - 6.11 Select Option 5 “MANUAL” from the Main Menu.
 - 6.12 Select Program 3.
 - 6.13 Select RINSE from the Manual Mode Status.
 - 6.14 Repeat Rinse a total of 3 times.
- 6.2 Drying Capillary with Air
 - 6.21 Place an empty sample vial in slot #1.
 - 6.22 Select Option 5 “MANUAL” from the Main Menu.
 - 6.23 Select Program 20.
 - 6.24 Select INJECT from the Manual Mode Status.