Western Regional Laboratory Drug Chemistry Procedure Manual Effective Date: June 19,1997

Name of Procedure:

Preliminary Tests Marquis Reagent

Suggested Uses:

The Marquis reagent consists of a solution of formaldehyde in concentrated sulfuric acid. Aromatic compounds that typically undergo electrophilic substitution will react with the Marquis reagent to produce colored intermediates. A positive response with the Marquis reagent is indicated by a significant color formation within 1-2 minutes.

Apparatus Needed to Perform Procedure Including Preparation of Reagent:

Fume hood

Gloves

Eye protection

Laboratory coat

Pipet with bulb

Graduated cylinder

50ml beaker

Glass stirring rod

Sulfuric acid (concentrated)

Formaldehyde solution (40%)

Funnel

Reagent bottle

Porcelain spot plate

Spatula

Formula for Preparing Reagent:

- 1. Measure out 10 milliliters of concentrated sulfuric acid in a beaker.
- 2. Add 8-10 drops of formaldehyde solution (40%) and stir.
- 3. Pour solution into a reagent bottle.
- 4. Properly label reagent bottle.

Formula for Preparing Reagent (continued):

Western Regional Laboratory Drug Chemistry Procedure Manual Effective Date: June 19,1997

Alternate Method

- 1. Pour 15-20 milliliters of concentrated sulfuric acid into a reagent bottle.
- 2. Add 0.2 0.3 gram of trioxane (trioxymethylene) and stir until completely dissolved.
- 3. Properly label reagent bottle.

Quality Control Check:

A quality control check of this reagent will be performed using a known standard of heroin and following the application procedure listed below.

Expiration Date of Chemical Reagent:

The Marquis reagent should be prepared every 30 days.

Application of Procedure on Evidence:

- 1. Place 1-2 drops of the reagent into a clean well on a spot plate.
- 2. With a spatula, add approximately 0.1 milligram of the unknown powder/tablet to the reagent in the spot plate.
- 3. Observe color produced.
- 4. Record results.

Safety Concerns:

Always wear eye protection, gloves and a laboratory coat when preparing this reagent.

Eye protection and a laboratory coat should be worn when using this reagent for color tests.

Sulfuric acid is a strong oxidizing agent and corrosive.

Western Regional Laboratory Drug Chemistry Procedure Manual Effective Date: June 19,1997

Literature References:

Wingrove, A.S., and Caret, R.L., <u>Organic Chemistry</u>, Harper and Row, New York, 1981, pp. 793-794.

Noller, C.R., <u>Chemistry of Organic Compounds</u>, W.B. Saunders Company, Philadelphia, 1965, p. 244.

Feigl, F., Spot Tests in Organic Analysis, 5th ed., Elsevier, New York, 1956, pp. 133-36.

Moffat, A. C., ed., <u>Clarke's Isolation and Identification of Drugs</u>, 2nd Ed., Pharmaceutical Press, London, 1986, p. 139.

Johns, S. H., "Spot Tests: A Color Chart Reference for Forensic Chemists," <u>Journal of Forensic Science</u>, July 1, 1979, pp. 631-649.

Butler, William P., <u>Methods of Analysis</u>, IRS Publication #341, December 1966, p. 136.