Western Regional Laboratory Drug Chemistry Procedure Manual Effective Date: June 19,1997

Name of Procedure:

Preliminary Tests Koppanyi Reagent

Suggested Uses:

The Koppanyi reagent consists of filter paper treated with (soaked in) cobalt acetate, methanol and glacial acetic acid and a solution of 5% isopropylamine in methanol. Barbiturates (5,5-disubstituted malonylureas) and other compounds will produce a positive color formation in less than 30 seconds.

Apparatus Needed to Perform Procedure Including Preparation of Reagent:

Fume hood Gloves Eye protection Laboratory coat Pipet with bulb Graduated cylinder 50ml beaker Glass stirring rod Filter paper Scissors Glacial acetic acid Cobalt acetate Methanol Wide mouth bottle with top Funnel Reagent bottle Porcelain spot plate

Spatula

Formula for Preparing Reagent:

- 1. Dissolve 0.1 gram of cobalt acetate in 100 milliliters of methanol.
- 2. Add 0.2 milliliter of glacial acetic acid.
- 3. Soak filter paper in the solution and allow to completely dry.

Western Regional Laboratory Drug Chemistry Procedure Manual Effective Date: June 19,1997

Formula for Preparing Reagent (continued):

- 4. Cut filter paper into approximately 1 inch squares.
- 5. Store filter paper in a wide mouth bottle with top.
- 6. Make a separate solution of 5% isopropylamine in methanol.
- 7. Pour into reagent bottle.
- 8. Properly label reagent bottles.

Quality Control Check:

A quality control check of this reagent will be performed using a known standard of a barbiturate and following the application procedure listed below.

Expiration Date of Chemical Reagent:

No expiration date. Reagents need to be properly contained in a sealed container and stored in a cool place.

Application of Procedure on Evidence:

- 1. Place a small amount (approximately 0.1 milligram) of sample on a piece of the Koppanyi paper.
- 2. Press the sample into the paper well with a spatula.
- 3. Place a drop of 5% isopropylamine solution on the edge of the Koppanyi paper and allow the drop to meet the sample.
- 4. Observe the color produced.
- 5. Record results.

Western Regional Laboratory Drug Chemistry Procedure Manual Effective Date: June 19,1997

Safety Concerns:

Always wear eye protection, gloves and a laboratory coat when preparing this reagent.

Eye protection and a laboratory coat should be worn when using this reagent for color tests.

Literature References:

Butler, William P., Methods of Analysis, IRS publication #341, December 1966, p.106.