Western Regional Laboratory
Drug Chemistry Procedure Manual
Effective Date: June 19,1997

Name of Procedure:

Thin-Layer Chromatography lodoplatinate Visualizing Reagent

Suggested Uses:

A visualizing reagent or detection reagent must be used in Thin-Layer Chromatography if the compound or compounds are not distinguishable by their own color. Opium alkaloids, coca alkaloids and most nitrogen containing compounds can be visualized using the iodoplatinate spray reagent.

<u>Apparatus Needed to Perform Procedure Including Preparation of Reagent:</u>

Fume hood

Graduated cylinder

Eye protection

Balance

Laboratory coat

Gloves

Spray bottle

Air compressor

Funnel

Spatula

Bottles

Tygon or rubber tubing

Platinic chloride (chloroplatinic acid)

Potassium iodide

Water

Formula for Preparing Reagent:

- 1. Dissolve 1 gram of platinic chloride (chloroplatinic acid) to 10 milliliters of water.
- 2. Add this solution to 10 grams of potassium iodide dissolved in 250 milliliters of water.
- 3. Dilute the solution from Step 1 to 500 milliliters with water.
- 4. Place desired amount in spray bottle.

Western Regional Laboratory Drug Chemistry Procedure Manual Effective Date: June 19,1997

Formula for Preparing Reagent (continued):

- 5. The remaining iodoplatinate solution can be stored in a reagent bottle in a cool place.
- 6. Properly label spray reagent.

Expiration Date of Chemical Reagent:

The iodoplatinate reagent can be used until depletion.

Application of Procedure on Evidence:

- 1. Place well-dried TLC plate in hood.
- 2. Activate hood.
- 3. Using the air compressor and spray bottle, apply a fine mist of the visualizing reagent to the TLC plate.
- 4. Apply the visualizing reagent until the spot corresponding to the known standard appears.
- 5. Compare the known standard and the compound in question for their size, shape, color and position on the TLC plate.
- 6. Record the results of your observation.

Safety Concerns:

Always wear eye protection, gloves, and a laboratory coat when preparing this reagent for use.

Eye protection and a laboratory coat should be worn when visualizing the TLC plate.

Western Regional Laboratory Drug Chemistry Procedure Manual Effective Date: June 19,1997

Literature References:

Randerath, Kurt, **Thin Layer Chromatography**, New York, Academic Press, 1968.

Moffat, A.C., <u>Clarke's Isolation and Identification of Drugs</u>, 2nd Ed., The Pharmaceutical Press, 1986.

Butler, William P., Methods of Analysis, IRS Publication #341, December 1966.