Western Regional Laboratory
Drug Chemistry Procedure Manual
Effective Date: June 19,1997

Name of Procedure:

Thin-Layer Chromatography p-DMAB Visualizing Reagent

Suggested Uses:

A visualizing reagent or detection reagent must be used in Thin-Layer Chromatography if the compound or compounds are not distinguishable by their own color. Suggested uses for the p-DMAB spray reagent are indole compounds, ergot alkaloids and other organic compounds.

Apparatus Needed to Perform Procedure Including Preparation of Reagent:

Fume hood

Graduated cylinder

Eye protection

Balance

Laboratory coat

Gloves

Spray bottle

Air compressor

Funnel

Spatula

Bottles

Tygon or rubber tubing

p-DMAB (para-dimethylaminobenzaldehyde)

Ethanol

Methanol

Concentrated hydrochloric acid

Formula for Preparing Reagent:

- 1. Weigh out 0.5 gram of p-DMAB (para-dimethylaminobenzaldehyde).
- 2. Add 5 milliliters of concentrated hydrochloric acid to 95 milliliters of methanol or ethanol.
- 3. Add the 0.5 gram of p-DMAB to the alcohol-hydrochloric acid mixture.

Formula for Preparing Reagent (continued):

Western Regional Laboratory Drug Chemistry Procedure Manual Effective Date: June 19,1997

- 4. Place this solution in the spray bottle.
- 5. Properly label spray reagent.

Expiration Date of Chemical Reagent:

The p-DMAB spray reagent is stable and can be used until depletion.

Application of Procedure on Evidence:

- 1. Place well-dried TLC plate in hood.
- 2. Activate hood.
- 3. Using the air compressor and spray bottle, apply a fine mist of the visualizing reagent to the TLC plate.
- 4. Apply the visualizing reagent until the spot corresponding to the known standard appears.
- 5. Compare the known standard and the compound in question for their size, shape, color and position on the TLC plate.
- 6. Record the results of your observation.

Safety Concerns:

Always wear eye protection, gloves, and a laboratory coat when preparing this reagent for use.

Eye protection and laboratory coat should be worn when visualizing the TLC plate.

Western Regional Laboratory Drug Chemistry Procedure Manual Effective Date: June 19,1997

Literature References:

Randerath, Kurt, **Thin Layer Chromatography**, New York, Academic Press, 1968.

Moffat, A.C., <u>Clarke's Isolation and Identification of Drugs</u>, 2nd Ed., The Pharmaceutical Press, 1986.

Butler, William P., Methods of Analysis, IRS Publication #341, December 1966.