Western Regional Laboratory Drug Chemistry Procedure Manual Effective Date: June 19,1997

Name of Procedure:

Extractions and Separations

Extraction of Amphoteric Compounds Including Morphine, Hydrocodone and other Alkaloids

Suggested Uses:

This is a general extraction procedure used to isolate and purify various amphoteric compounds.

Apparatus Needed to Perform Procedure Including Preparation of Reagent:

Fume hood

Gloves

Eye protection

Laboratory coat

Pipet with bulb

Graduated cylinder

Glass stirring rod

pH Test paper

Separatory funnel

Funnel

Sodium Bicarbonate

Chloroform

Isopropanol

Ethyl Ether

5% Hydrochloric Acid Reagent

Formula for Preparing Reagent:

5% Hydrochloric Acid Reagent

- 1. Measure 95 milliliters of water in a 100 milliliter graduated cylinder.
- 2. Bring to total volume (100ml) with concentrated hydrochloric acid.
- 3. Pour into a reagent bottle.
- 4. Properly label reagent bottle.

Formula for Preparing Reagent (continued):

Western Regional Laboratory Drug Chemistry Procedure Manual Effective Date: June 19,1997

3:1 Chloroform/Isopropanol Reagent

- 1. Measure out 60 milliliters of chloroform.
- 2. Measure out 20 milliliters of isopropanol.
- 3. Combine the chloroform and isopropanol.
- 4. Pour solution into a reagent bottle.
- 5. Properly label reagent bottle.

Expiration Date of Reagent:

Reagent may be used until depletion if stored In an airtight reagent bottle.

Application of Procedure on Evidence:

- 1. Dissolve sample in 5% hydrochloric acid reagent.
- 2. Extract with ethyl ether and discard ether washings.
- 3. Add sodium bicarbonate until a pH of 8 is obtained. Check with pH test paper.
- 4. Extract the basic solution with the 3:1 chloroform/isopropanol reagent.
- 5. The solvent extract may be dried using magnesium sulfate of sodium sulfate.
- 6. Evaporate the solvent using dry nitrogen and heat to obtain the free base form of the compound.

Safety Concerns:

Ethyl Ether and Isopropanol are flammable solvents. Chloroform should be used in a fume hood.

Literature References:

Western Regional Laboratory Drug Chemistry Procedure Manual Effective Date: June 19,1997

Shriner, R.L., Fuson, R. C., Curtin, D. Y., <u>The Systematic Identification of Organic Compounds</u>, 5th Ed., Wiley & Sons Inc., New York 1964, pp. 67-85, 88-106.

Butler, William P., Methods of Analysis, IRS Publication #341, December 1966, p. 64.