

**Name of Procedure:**

Extractions and Separations  
Dry Solvent Extraction of Drugs Utilizing Chloroform/Ammonia

**Suggested Uses:**

This procedure uses a dry extraction of chloroform saturated with ammonia to remove hydromorphone, morphine, diazepam, lorazepam, flurazepam, phentermine, chlordiazepoxide, cocaine, pentazocine, methaqualone, and benzodiazepines from pharmaceutical preparations and clandestine mixtures.

**Apparatus Needed to Perform Procedure Including Preparation of Reagent:**

Fume hood  
Gloves  
Eye protection  
Laboratory coat  
Ammonium Hydroxide  
Chloroform  
Small beaker  
Filter paper  
Pipets, glass, disposable  
Pipet bulb  
Heat source  
Reagent bottle

**Formula for Preparing Reagent:**

1. Mix approximately one part ammonium hydroxide to ten parts chloroform and place in a reagent bottle.
2. Shake reagent bottle and allow layers to separate.
3. Properly label bottle.

**Expiration Date of Chemical Reagent:**

The reagent can be used until depleted if stored in an airtight reagent bottle.

**Application of Procedure on Evidence:**

1. Place 20-50 milligrams of sample in filter paper over a small beaker.
2. Wash sample with ethyl ether and discard washings.
3. Dry the sample and then wash with several small portions of chloroform/ammonia reagent.
4. Evaporate solvent over moderate heat in a fume hood.

**Safety Concerns:**

Ammonia is a strong base and is caustic. Ammonia and chloroform should be used in a well-ventilated area or under a fume hood. Ethyl ether is very flammable.

**Literature References:**

Adair, A., Noggle, F. Jr., Odom, M., Rhodes, M., "The ANOR (Alternate Non-aqueous Organic Ratio Extraction Procedure)", **MICROGRAM**, Vol. XVI., No. 1, Jan. 1983, pp. 220-224.